For better dentistry

**DENTSPLY**

# AquaCem<sup>®</sup>

**Glass Ionomer Luting Cement**

**Glasionomer Befestigungszement**

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CAUTION: For dental use only.

USA: Rx only.

Content.....	Page
<b>1 Product description.....</b>	<b>2</b>
<b>2 Safety notes .....</b>	<b>3</b>
<b>3 Step-by-step instructions .....</b>	<b>4</b>
<b>4 Lot number and expiration date.....</b>	<b>5</b>

### 1 Product description

**AquaCem<sup>®</sup>** is a light yellow, translucent glass ionomer luting material consisting of a blend of alumino-silicate glass and polyacrylic acid. The powder is mixed with distilled water to produce a luting material which adheres to dentin and enamel producing tightly sealed cementations.

**AquaCem<sup>®</sup>** complies with the requirements of ISO 9917 for polyalkenoate luting cements.

#### 1.1 Delivery forms

- A powder/liquid version for hand mixing.

#### 1.2 Composition

- Calcium-sodium-fluoro-phosphoroaluminium-silicate
- Polyacrylic acid
- Tartaric acid
- Iron oxide pigments

#### 1.3 Indications for use

- Cementation of crowns, inlays, bridges, and orthodontic bands.

#### 1.4 Contraindications

- Direct or indirect application to dental pulp (direct or indirect pulp-capping).
- Cementation of extensive bridgework when the time needed for application of the cement and positioning of the bridge exceeds the working time (2 minutes 30 seconds) of AquaCem<sup>®</sup>.

## 2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

### Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

### 2.1 Warnings

Components and mixed material may be irritating to skin, eyes and oral mucosa.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes rinse with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, remove material from the tissues. Flush mucosa with plenty of water and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

### 2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- To prevent bottles, powder scoop and mixing pad from exposure to spatter or spray of body fluids or contaminated hands it is recommended that the bottles, powder scoop and mixing pad are handled offside the dental unit with clean/disinfected gloves.
- Tightly close bottle immediately after use.
- The desiccant crystals integrated in the cap and/or drying capsule protect the powder from absorbing moisture. The color of the desiccant crystals should be orange. The desiccant crystals are no longer effective after the color has turned from orange to colorless.
- Increased temperature and/or humidity reduce available working time and accelerate setting time. Material will set faster in the mouth than at the bench. Do not use bench cure rates to judge mouth cure rates.
- The product should be mixed in proportions of one level scoop powder to two drops water<sup>1</sup> for optimum performance. The mixed material should be homogeneous and streak free prior to application. Variations may affect the strength and durability of the material.
- Use in a well-ventilated area. Avoid inhaling dust particles from powder.

### 2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation. Reddish rashes may be seen on the skin.
- Contact with mucous membranes: Inflammation (see Warnings).

### 2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store with original (not dropper cap) caps tightly closed in a well-ventilated place at temperatures between 10 °C and 24 °C (50 °F and 75 °F).
- Allow material to reach room temperature prior to use.
- Protect from moisture.
- Do not use after expiration date.

## 3 Step-by-step instructions

### 3.1 Cavity preparation

1. Isolate the field of operation in preferred manner. Rubber dam is recommended.
2. Cavity preparation must provide for the mechanical retention of the product.
3. If caries is present, completely excavate using low-speed and/or hand instrumentation.
4. Place matrix and wedge where indicated.
5. Gently dry preparation with cotton pellet. Avoid desiccation.

### 3.2 Pulp protection

1. In deep cavities cover the dentin close to the pulp (less than 1 mm) with a hard-setting calcium hydroxide liner (Dycal® calcium hydroxide liner, see complete Instructions for Use).

### 3.3 Dosage and mixing

1. Invert powder bottle to fluff the powder.
2. Overfill the scoop with powder and strike off the excess with the bottle insert leaving a flat surface level with the edge of the scoop.

### NOTICE

#### Premature setting of material

Reduced working time.

1. Keep powder and liquid separated until ready for mixing.
3. Place the powder onto the pad provided or on a glass slab<sup>2</sup>. Special spatulas are not required.
4. Divide the powder into two halves.
5. Fill water dispenser with distilled or de-ionized water.
6. Dispense two drops of water by holding the bottle vertically and squeezing gently.
7. Add one half of the powder to the water. As soon as the first half is wetted, incorporate the second, then mix thoroughly. The measures are intended as a guide; the mix should have a smooth creamy consistency. It is essential that the powder is incorporated very quickly into the water to ensure the optimum mix ratio. For all normal uses the mixing time should not exceed 15 seconds<sup>2</sup>.

### 3.4 Working and setting time<sup>3</sup>

- From start of mixing, the remaining working time is approximately 2 minutes 30 seconds.
- After mixing, the setting time is 3 minutes 30 seconds to 5 minutes.

### 3.5 Application

1. Apply cement onto the fitting surface of the restoration/orthodontic band and/or the prepared tooth immediately upon completion of the mix whilst appearance is glossy.
2. Position the restoration/orthodontic band and apply adequate pressure to seat it. Avoid contamination with water/saliva during application and setting. It should, however, be noted that the cement remains subject to dehydration if left isolated for longer than 10 minutes from start of mixing. This means that water balance must be maintained by releasing the cement to the oral environment within this time range.

### 3.6 Finishing

1. Excess material may be removed at once or after completion of setting time (see 3.4 Working and setting time). Do not touch the cement when it is in the intermediate, rubbery condition. For removal of surplus material, appropriate hand instruments, like a dental probe, and dental floss are used.

### 3.7 Cleaning instruments

Glass ionomers adhere to instruments and glass slabs.

1. Immediately after use, either clean instruments with water or drop in disinfectant or water to facilitate subsequent cleaning.

## 4 Lot number ( ) and expiration date ( )

1. Do not use after expiration date. ISO standard is used: "YYYY-MM"
2. The following numbers should be quoted in all correspondence:
  - Reorder Number
  - Lot number
  - Expiration date

© DENTSPLY DeTREY 2013-04-17

[These Instructions for Use are based on Master Version 02]

---

<sup>1</sup> For tests according to ISO 9917 the recommended powder : liquid ratio is 3.3 g : 1 g by weight.  
<sup>2</sup> Mixing and application of the material should be done at normal room temperature. Higher temperatures shorten the working and setting times of the cement. Under warm conditions or if a long working time is required, it is therefore recommended that the glass slab be cooled. Testing to ISO 9917 is carried out at 23 °C ± 1 °C and a relative humidity of 50 ± 10%.  
<sup>3</sup> Net setting time determined according to ISO 9917.

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**DENTSPLY DeTrey GmbH**
**DENTSPLY Limited
**DENTSPLY Italia S.r.l.