Esthetic Resin Cement

DIRECTIONS FOR USE – ENGLISH

For dental use only.
USA: Rx only.

1. PRODUCT DESCRIPTION

Calibra® Esthetic Resin Cement is a visible light cured, dual cured or self cured high strength resin cement compatible with numerous dentin/enamel adhesive systems. It is to be used following the application of a suitable dentin/enamel adhesive to adhesively bond and lute indirect restorations to tooth structure.

1.1 Delivery forms*

Calibra® Cement and Try-In pastes are available in:

- small syringes for easy dispensing
- a convenient dual-barreled syringe
- 5 shades: translucent, light, medium, dark, and opaque

*Some delivery forms may not be available in all countries

1.2 Composition

- Calibra® Esthetic Resin Cement Base: Dimethacrylate Resins; Camphorquinone (CQ) Photoinitiator; Stabilizers; Glass Fillers; Fumed silica; Titanium Dioxide; Pigments
- Calibra® Esthetic Resin Cement Catalyst: Dimethacrylate Resins; Peroxide Catalyst; Stabilizers; Glass Fillers; Fumed Silica
- Calibra® Esthetic Resin Cement Try-In Paste: Glycerine; Fumed Silica; Titanium Dioxide; Pigments
- Calibra® Silane Coupling Agent: Acetone; Ethyl Alcohol; Organo Silane

1.3 Indications

Adhesive Cementation of:

1. Ceramic, porcelain, composite inlays/onlays, veneers and crowns.
2. All metal crowns, bridges, inlays/onlays including precious, semi-precious and non-precious metals.
3. PFM (porcelain fused to metal) crowns and bridges.
4. Prefabricated and cast posts.
5. Resin-bonded retainer bridges (Maryland bridges).

1.4 Contraindications

1. Calibra® Cement is contraindicated for use with patients who have a history of severe allergic reaction to methacrylate resins or any of the components.

1.5 Compatible adhesives

Calibra® Cement is used following application of a suitable dentin/enamel adhesive and is chemically compatible with conventional methacrylate-based dentin/enamel adhesives including DENTSPLY adhesives designed for use with visible light cured composite restoratives for Calibra® Cement light cured applications and DENTSPLY adhesives designed for use with dual-cured resin based materials for Calibra® Cement dual cured applications. See complete directions for use of selected adhesive.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other sections of these directions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Calibra® Cement contains polymerizable methacrylate monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.
 - **Avoid Eye Contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
 - **Avoid Skin Contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs immediately remove material with cotton and alcohol and wash well thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - **Avoid Contact with Oral Soft Tissues/Mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
2. Calibra® Silane Coupling Agent contains acetone. Do not breathe vapor.

2.2 Precautions

1. This product is intended to be used only as specifically outlined in the Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and sole responsibility of the practitioner.

2. Insufficient data exist to support the use of Calibra® Esthetic Resin Cement as a core build up, filling material, base or liner.
3. Insufficient data exist to support the use of Xeno® Self-Etch Adhesives (all versions) with Calibra® Esthetic Resin Cement for veneer cementation light cured.
4. Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
5. Contact with saliva, blood and/or some astringent solutions during adhesive procedures may cause failure of the restoration. Use of rubber dam or adequate isolation is recommended.
6. Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
7. Increased temperature may accelerate the setting reaction. After placing Calibra® Cement in contact with tooth structure, e.g., within endodontic post space or in inlay/onlay preparations, immediately seat restoration. Any delay may allow polymerization to begin, which may prevent complete seating of the restoration.
8. Exercise caution when placing multiple units not to exceed available working time. Load and insure complete seating of each unit prior to loading successive units. The use of a fresh mix of cement is recommended if available working time is approached (see Interactions).
9. Calibra® Cement, when used with Catalyst for self-cure or dual cure should be mixed in equal volumes for optimum performance. Variations from this 1:1 ratio may shorten available working time and compromise final physical properties.
10. Containers should be tightly closed immediately after use. Syringe should be tightly closed by replacing the original cap immediately after use.
11. Calibra® Cement should extrude easily. DO NOT USE EXCESSIVE FORCE. Excessive pressure may result in unanticipated extrusion of the material or cause syringe rupture.
12. The Calibra® Cement is a light-cured material. Proceed immediately once materials have been placed on mixing pad or protect from ambient light.
13. Use only in well ventilated areas.
14. Flammable: Calibra® Silane Coupling Agent contains acetone. Keep away from sources of ignition.

Interactions

- Eugenol containing materials should not be used in conjunction with this product because they may interfere with hardening and cause softening of the polymeric components of the material.
- Contact with some astringent solutions may interfere with hardening of the polymeric components of the material.
- As with any dual cure resin cement system, the use of a dual cure adhesive system can shorten working time. This effect should be investigated in the laboratory prior to clinical use.
- Variable in-vitro data exist regarding use of Calibra® Esthetic Resin Cement in dual cured or self-cured (limited or no light curing) applications in conjunction with some light-cured-only adhesives. Chemical/Product incompatibility may adversely affect product efficacy, leading to premature restoration failure.

2.3. Adverse reactions

Product may irritate the eyes and skin. **Eye contact:** irritation and possible corneal damage. **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin. **Mucous membranes:** inflammation (see Warnings).

2.4. Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Keep out of direct sunlight and store in a well ventilated place at temperatures between 2°-24°C/36°-75°F. Allow material to reach room temperature prior to use. Protect from moisture. Do not freeze. Do not use after expiration date.

3. STEP-BY-STEP INSTRUCTIONS

The following preparatory steps apply to cementation of all types of indirect restorations:

- I. Following removal of the temporary restoration and any remaining temporary cement, clean enamel and dentin as directed by the adhesive manufacturer's directions.
- II. Check the fit and esthetics of the restoration. **Technique Tip:** Occlusal adjustment of inlays/onlays/veneers is best accomplished after final cementation.

TRY-IN PASTE (OPTIONAL)

1. Dispense the appropriate shade of try-in paste from the syringe onto a clean mixing pad or glass slab. Load paste onto internal surfaces of the restoration and gently seat onto preparation. Clean excess with a cotton pellet and/or blunt explorer. Shades may be blended to achieve optimum esthetics. **Important Technique Tip:** Factors such as restoration shape, texture, thickness and shade(s) of the underlying and surrounding tooth structure contribute to final perceived shade. Try-in paste is a guide for cement shade range selection only. **NOTE:** The try-in paste will not polymerize, thus offers unlimited work time.
2. Once fit and esthetics are verified, thoroughly rinse try-in paste from restoration and preparation surfaces using water.

TREATMENT OF THE RESTORATION

Metal Restorations

Internal surfaces of restorations should be clean and dry prior to cementation. Internal surface microetching (sandblasting with 50µ alumina) of metal surfaces of the restoration is recommended.

Technique Tips For Maryland Bridges: The use of 180° wrap of wings, rest seats, parallelism and slots or grooves is necessary for appropriate retention. The metal wings of the Maryland Bridge should be perforated, electrolytically etched, laboratory chemically etched, or mechanically sandblasted with 50µ alumina.

Ceramic/Composite Restorations

Follow the dental laboratory or restoration manufacturer's instructions for pre-treatment, if required. Restorations designed to be silanated or if the internal silanated surface has been disturbed during try-in, apply Calibra® Silane Coupling Agent as outlined below:

1. Prior to applying the silane, clean the internal surface of the restoration with Caulk® 34% Tooth Conditioner Gel (34% phosphoric acid) (see complete Directions for Use). Apply for 30 seconds and rinse thoroughly with water for 20 seconds. Air-dry the internal surface of the restoration with clean, dry air from a dental syringe.

Danger of injury due to excessive force

- Apply slow and steady pressure on the syringe plunger
- Do not use excessive force – rupture may result.

2. Attach supplied needle tip to end of the Calibra® Silane Coupler syringe. Gently apply pressure to syringe plunger. DO NOT USE EXCESSIVE FORCE. Silane agent should express easily, one drop at a time. Directly apply the silane to the clean internal surface of the restoration. Air-dry the silanated surface of the restoration with clean, dry air from a dental syringe or allow to air dry. **NOTE:** Avoid application of the silane to the external surface of the restoration. Contact with the external surface will cause bonding of the cement to the external surface, complicating clean-up and necessitating removal of the glazed portion of the external surface of the restoration.
3. If silanated surface becomes contaminated, clean surface with alcohol or acetone, allow to air dry and repeat application as outlined above.

TOOTH CONDITIONING/DENTIN PRETREATMENT/ADHESIVE APPLICATION

1. Following restoration try-in, rinse preparation thoroughly with water spray and air dry.
2. Proceed with adhesive application as directed in the complete directions for use.

Delayed or Reduced Adhesion due to adhesive incompatibility

- Select compatible light cure or dual cure adhesive
- Follow selected adhesive directions for use

3. Always refer to and follow adhesive manufacturer's complete directions for use for selection and application of the appropriate light cure or dual cure adhesive. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately following adhesive application to cementation with Calibra® Cement, light cured base or dual or self cured with catalyst as indicated in the appropriate section below.

1. Veneer cementation – light-cured

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment/Adhesive Application. Adhesives compatible with light-cured, CQ initiated methacrylates are recommended. If veneer restorations exceed 1.0mm thickness, are heavily opaqued, or otherwise impede the transmission of light, operator is urged to follow section 2.0, Dual-Cured cementation technique. Apply adhesive to the internal bonding surface of the restoration if directed by adhesive manufacturer's directions for use.

1.1 Cementation technique

- 1.1.1 Dispense the desired shade of Calibra® Esthetic Resin Cement base paste from the syringe directly onto the veneer. Protect cement from exposure to light.
- 1.1.2 Seat the loaded veneer in place. Remove excess from the gingival margin with a blunt instrument. Briefly (10 seconds or less) light cure the gingival portion only to tack restoration in place. Remove any excess from proximal and lingual margins. **Technique Tip:** Placement of mylar strips between preparation and adjacent teeth prior to seating veneer aids in isolation and in excess cement clean up. After removing gingival excess and tacking veneer in place, remove interproximal excess cement by pulling mylar strip towards the facial, i.e., from tooth to restoration surface.

Inadequate polymerization due to insufficient curing

- Check compatibility of curing light
- Check curing cycle
- Check curing output before each procedure

- 1.1.3 Using a visible light, curing unit designed to cure CQ initiated methacrylates (spectral output including 470nm), with a minimum output of 550mW/cm², light cure 20 seconds each from the buccal, lingual, and interproximal aspects.
- 1.1.4 Following the light cure check and adjust occlusion as necessary. Proceed to Finishing and Polishing section.

2. Crown/bridge/inlay/onlay cementation-dual cured

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment/Adhesive Application. Adhesives (plus activators, if applicable) compatible with both light-cured, CQ initiated methacrylates and peroxide-amine initiated self-cured methacrylates are recommended. Apply adhesive/activator (if applicable) to the internal bonding surface of the restoration if directed by adhesive manufacturer's directions for use.

2.1 Cementation technique

- 2.1.1 Dispense the desired shade of Calibra® Esthetic Resin Cement base paste from the syringe onto a clean mixing pad. Dispense an equal amount of the desired viscosity of catalyst paste (Regular or High Viscosity). Mix the cement for 20-30 seconds.

Technique Tip: As with any dual cure resin cement system, the use of a dual cure adhesive system can shorten working time (See Interactions). Calibra® Esthetic Resin Cement base may be mixed with High Viscosity catalyst yielding a High Viscosity, "stay-put" thixotropic dual cured cement suitable for inlays, shallow onlays and veneers. Mixing Calibra® Esthetic Resin Cement base with Regular Viscosity catalyst yields an easier flowing dual cured cement suitable for large, multi-surface and full coverage restorations with higher hydraulic seating pressure.

Danger of injury due to excessive force

- Apply slow and steady pressure on the syringe
- Do not use excessive force – dual barrel syringe rupture may result

Dual Barreled syringe dispensed

1. Remove syringe cap. Dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess so base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
 2. Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mixing tip with the v-shape notch on the syringe flange. Turn colored mixing tip cap 90 degrees in a clockwise direction to lock in place on syringe.
 3. Gently depress syringe plungers to begin the flow of material. DO NOT USE EXCESSIVE FORCE. If force is encountered, remove syringe from operating field, remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.
- 2.1.2 Apply a uniform layer of cement on the entire internal surface of the restoration. For inlays/onlays, it may be helpful to apply a thin layer of cement to the internal portions of the tooth preparation to avoid any porosity or voids. **Technique Tip:** At room temperature, Calibra® Esthetic Resin Cement, in either viscosity, offers a minimum work time of 2 min. 30 sec. on the mixing pad when protected from ambient light. If loaded into restoration lined with adhesive/activator mixture, room temperature work time (restoration seating time) will be less (see Precautions).
 - 2.1.3 Seat the restoration with gradual pressure. A gentle rocking or vibratory motion may be helpful to insure optimal seating.
 - 2.1.4 Remove gross excess from marginal areas. Use an instrument such as a blunted explorer, periodontal probe or a clean, dry brush tip. Restoration should not be moved or torqued during removal of gross excess cement. Special attention should be paid to interproximal areas using floss to remove excess cement. **Technique Tip:** A 10 second light "pre-cure" of excess cement at the margins will cause cement to "gel", allowing easy cleanup. **Technique Tip:** Apply moderate and consistent pressure to the restoration throughout the self-cure set time of approximately 6 minutes from the beginning of mixing.
 - 2.1.5 Once stabilized, light cure, using a visible light, curing unit designed to cure CQ initiated methacrylates (spectral output including 470nm), with a minimum output of 550mW/cm². Cure all marginal areas of the restoration for 20 seconds from each direction – buccal, lingual and occlusal.
 - 2.1.6 Following the self-cure set check and adjust occlusion as necessary. Proceed to Finishing and Polishing, section.

3. Endodontic posts

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment/Adhesive Application. Adhesives (plus activators, if applicable) compatible with both light-cured, CQ initiated methacrylates and peroxide-amine initiated self-cured methacrylates are recommended. Apply adhesive/activator (if applicable) to the bonding surface of the post if directed by adhesive manufacturer's directions for use.

3.1 Cementation technique

- 3.1.1 Dispense the desired shade of Calibra® Esthetic Resin Cement base paste from the syringe onto a clean mixing pad. A light base shade e.g., Translucent will allow maximum light transmission. Dispense an equal amount of Regular Viscosity catalyst paste. Mix the cement for 20-30 seconds. **Technique Tip:** As with any dual cure resin cement system, the use of a dual cure adhesive system can shorten working time. (See Interactions)

Danger of injury due to excessive force

- Apply slow and steady pressure on the syringe
- Do not use excessive force – dual barrel syringe rupture may result

Dual Barreled syringe dispensed

1. Remove syringe cap. Dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess so base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
 2. Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mixing tip with the v-shape notch on the syringe flange. Turn colored mixing tip cap 90 degrees in a clockwise direction to lock in place on syringe.
 3. Gently depress syringe plungers to begin the flow of material. DO NOT USE EXCESSIVE FORCE. If force is encountered, remove syringe from operating field, remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.
- 3.1.2 Spread mixed Calibra® Esthetic Resin Cement components on surface of post and/or into the post preparation with a syringe tip, Lentulo Spiral or metal file.
 - 3.1.3 Seat post immediately. Clean up excess with appropriate instruments. A 10 second light exposure "pre-cure" of excess cement at the margins will cause cement to "gel", allowing easy cleanup.

- 3.1.4 Stabilize post until cement sets. Apply moderate and consistent pressure to the restoration throughout the self-cure set time of approximately 6 minutes from the beginning of mixing.
- 3.1.5 Once stabilized, light cure, using a visible light, curing unit designed to cure CQ initiated methacrylates (spectral output including 470nm), with a minimum output of 550mW/cm². Cure all accessible areas of the post for 20 seconds.
- 3.1.6 Proceed with core build-up and/or preparation.

FINISHING AND POLISHING

1. Removal of resin flash is best accomplished with the Enhance® Finishing System of points, cups and discs. The Enhance® System will remove flash and finish restoration margins without removal of or trauma to the enamel.
2. Polish final restoration using PoGo® Polishing system or Prisma®•Gloss™ Polishing Paste and Prisma®•Gloss™ Extra Fine Polishing Pastes. (See complete Directions for Use supplied with polishing product chosen)

4. HYGIENE

Cross-contamination

- Do not reuse single use products. Dispose in accordance with local regulations.
- Reprocess reusable products as described below.

4.1 Cleaning and disinfection

- For dual-barreled syringe, remove used mixing tip and discard appropriately. Replace original syringe cap prior to storing. The dual-barreled syringe may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.
- To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid package contamination.
- Disinfect contaminated syringes with a water-based hospital-level disinfection solution according to national/local regulations.
- Repeated disinfection may damage label. **NOTE:** Destruction of the label by vigorous wiping. Wipe syringe gently.
- The supplied mixing spatula may be cleaned by scrubbing with hot water and soap or detergent. Do not autoclave spatula. Disinfect as outlined below.

4.2 Disinfection of re-usable mixing spatula

Disinfect spatula with a hospital-level, tuberculocidal disinfectant solution according to national/local regulations. Iodophors, sodium hypochlorite (5.25%), chlorine dioxide and dual or synergized quaternary ammoniums are approved disinfectants. Some phenolic-based agents and iodophor-based products may cause surface staining. The disinfectant manufacturer's directions should be followed properly for optimum results. Water-based disinfectant solutions are preferred.

5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. ISO standard is used: "YYYY/MM"
2. The following numbers should be quoted in all correspondence:
 - Reorder Number
 - Lot number
 - Expiration date

Manufactured by:
DENTSPLY Caulk
38 West Clarke Avenue
Milford, DE 19963 USA
Tel.: 1-302-422-4511
www.dentsply.com

DENTSPLY DeTREY GmbH
De-Trey-Str. 1
78467 Konstanz
Germany
Tel.: 49-7531-583-0
www.dentsply.de

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DENTSPLY DeTrey GmbH