For better dentistry

DENTSPLY

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CAUTION

Cross-contamination

After intraoral use do not clean, disinfect or reuse Compules® tips.

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DENTSPLY

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Standard Posterior Restorative

Dyract® Posterior is a compomer-based restorative material specially designed for posterior use.

Dyract Posterior exhibits advanced performance regarding abrasion resistance and physical properties which are of special importance for occlusal stress-bearing restorations.

Dyract Posterior allows a simplified and fast application technique and combines fluoride-releasing glass-ionomer chemistry with the strength and aesthetics of a composite.

Dyract Posterior, pre-dosed in Compules® Tips for direct intra-oral application, is shaded white-opaque for easy differentiation between tooth substance and restoration.

Dyract Posterior is used following application of Xeno®III Single Step Self-Etching Adhesive or Prime&Bond® NT, a universal self-priming dental adhesive designed to bond the restorative to enamel and dentine.

Caution: For dental use only.

COMPOSITION

Urethan dimethacrylate (UDMA)
Carboxylic acid modified dimethacrylate (TCB Resin)
Alkanoyl-poly-methacrylate
Strontium-fluoro-silicate glass
Strontium fluoride
Photo initiators
Butyl hydroxy toluene

INDICATIONS

Dyract Posterior is indicated for Class I and II cavities when aesthetic considerations are not given priority. Cavity width must be less than $2/3$ of the intercuspal distance.

CONTRAINdicATIONS

- Direct or indirect pulp capping.
- Core build-up for full ceramic crowns.

- Use in patients with a known allergy to dimethacrylate resins or any other of the components.
- When contamination with saliva, blood, etc. cannot be avoided.

WARNINGS

1. Dyract Posterior contains methacrylates that may be irritating to the eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

2. Dyract Posterior may cause sensitisation by contact with skin or mucous membranes in susceptible persons. After accidental contact, wash immediately with plenty of soap and water or rinse with plenty of water.

If sensitisation has occurred, discontinue use.

PRECAUTIONS

1. Use a gentle, even motion when exerting pressure on the Compules tip applicator gun. Use of excessive force or sudden movement could create a potential hazard.
2. Keep out of sunlight. To be stored at temperatures between 10 °C and 24 °C.

3. Humidity can adversely affect the properties of unsealed Compules tips. Therefore keep Compules tips sealed in their blister pack until use.

INTERACTIONS WITH DENTAL MATERIALS

Eugenol containing dental materials should not be used in conjunction with this product because they may inhibit the curing of monomeric components.

ADVERSE REACTIONS

In isolated cases, reversible inflammatory changes of the oral mucosa have been reported after accidental contact with acrylate monomers.

STEP-BY-STEP INSTRUCTIONS

1. Cavity Preparation

In all classes of cavities, this may be kept to the minimum required for caries removal.

2. Cleaning

Cavity cleanliness is paramount for the development of adhesion.

In cases where no cavity preparation has been made, clean the tooth surface with a rubber cup and pumice or a prophy-paste like Nupro[®]. Preparing a fresh surface with a finishing bur will significantly increase bond strength to enamel.

Wash surface thoroughly with air/water spray.

Remove rinsing water by blowing gently with an air syringe or blot-dry with a cotton pellet.

Do not desiccate the dentine structure.

3. Pulp Protection

For direct or indirect pulp-capping protect the dentine close to the pulp (< 1 mm) with a hard-setting calcium hydroxide liner (e.g. Dycal[®]), leaving the remaining cavity surface free for bonding with the adhesive.

4. Conditioning and application of adhesive

Prior to the application of Dyract Posterior, the cavity has to be treated with Xeno III – Single Step Self-Etching Dental Adhesive. Alternatively, Prime&Bond NT – Nano-Technology Dental Adhesive – may be used. For application please follow the directions for use of the respective adhesive.

5. Placement of Dyract Posterior

Insert Compules Tip into the notched opening of the applicator gun barrel.

Dispense Dyract Posterior directly into the cavity preparation. In deep cavities, incremental placement and curing (in 3 mm layers or less) is recommended to minimise polymerisation shrinkage.

6. Curing

Cure each increment (3 mm or less) separately with a VLC dental polymerisation unit for 20 seconds¹. The tip of the light guide should be held as close as possible to the restoration during curing.

¹ Check curing light for minimum curing output of at least 500 mW/cm².

Important: Be sure to expose each area of the entire restoration to the curing light. Additionally, the restoration should be cured through lingual or buccal enamel walls.

7. Finishing

Begin finishing immediately after curing. Gross excess material may be removed with fluted finishing burs or diamonds. Finishing is best achieved by using Enhance™

Finishing and Polishing Discs and interproximal finishing and polishing strips. A high final lustre can be obtained by applying Prisma® Gloss™ and Prisma Gloss Extrafine Polishing Pastes.

MAINTENANCE OF APPLICATOR GUN

The applicator gun is sterilisable by autoclave or cold sterilisation solution following the manufacturers' instructions.

It is recommended that the applicator gun be disassembled for assured sterilisation. Partially close the applicator gun and place thumb under the rear portion of the hinge. Push upward and lift hinge separating the applicator gun, exposing the plunger. Remove residual compomer with a soft paper tissue and a suitable solvent (70 % alcohol).

To reassemble, insert plunger into applicator gun barrel, press components together and snap hinge mechanism in place.

STORAGE

Keep out of sunlight.

To be stored at temperatures between 10 °C and 24 °C. Humidity can adversely affect the properties of unsealed Compules tips. Therefore keep Compules tips sealed in their blister pack until use. Under normal ambient conditions, unsealed Compules tips stay usable for about 4 weeks.

Inadequate storage conditions will shorten the shelf life and may lead to a malfunction of the product.

BATCH NUMBER LOT AND EXPIRY DATE

The batch number should be quoted in all correspondence which requires identification of the product. Do not use after expiry date.

If you have any questions, please contact:

Manufacturer:
DENTSPLY DeTREY GmbH
De-Trey-Str. 1
78467 Konstanz
GERMANY
Phone +49 (0) 75 31 5 83-0

Distributor:
DENTSPLY Limited
Building 1
Aviator Park, Station Road
Addlestone, KT15 2PG
GREAT BRITAIN
Phone (0 19 32) 85 34 22

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© DENTSPLY DeTREY 2003-10-07