For better dentistry

DENTSPLY

Dyract[®] Seal

Compomer Pit & Fissure Sealant

Kompomer-Fissurenversiegler

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CAUTION: For dental use only.

USA: Rx only

Content	Page
1 Product description	2
2 Safety notes	3
3 Step-by-step instructions	5
4 Hygiene	7
5 Lot number and expiration date	7

1 Product description

Dyract[®] Seal fissure sealant is a light-curing dental pit & fissure sealant containing releasable fluoride. The fluoride offers an additional protection of the sealed pits and fissures.

Dyract[®] Seal fissure sealant is available in a clear and a white-opaque shade.

Dyract[®] Seal fissure sealant – when used in combination with NRC™ Non-Rinse Conditioner – offers a "water free" application technique, which is of particular advantage with treatment of small children.

Dyract[®] Seal fissure sealant – when used in combination with Prime&Bond[®] NT – allows a complete obliteration even of deep and narrow fissures (Total Seal™ Technique).

1.1 Delivery forms

- Small syringes containing 1.5 g (1 ml) fissure sealant to be used in combination with disposable steel canulas (applicator needles) for direct intra-oral application.

1.2 Composition

- Strontium-alumino-fluoro-phosphor-silicate glass
- Highly dispersed silicon dioxide
- Ammonium salt of phosphoric acid modified methacrylate resin
- Carboxylic acid modified methacrylate resin
- Diethyleneglycol dimethacrylate (DGDMA)
- Camphorquinone
- Ethyl-4-dimethylaminobenzoate
- Butylated hydroxy toluene (BHT)
- Titanium dioxide (white version only)

1.3 Indications

- Preventive sealing of pits and fissures in the primary and secondary dentition.

1.4 Contraindications

- Use with patients who have a history of severe allergic reaction to methacrylate resins.
- Not to be used in case of caries.

1.5 Compatible adhesives

The material is chemically compatible with methacrylate-based dentin/enamel adhesives including DENTSPLY adhesives designed for use with visible light cured composite restoratives (see complete Instructions for Use of selected adhesive).

For the complete obliteration of deep and narrow fissures (Total Seal™ Technique) the use of Prime&Bond® NT adhesive is specially recommended.

2 Safety notes

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

- This is the safety alert symbol. It is used to alert you to potential personal injury hazards.
- Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

The material contains methacrylates and polymerizable monomers which may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after fissure sealing is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
Remove soft (oxygen-inhibited) surface layer after light curing with cotton pellets or cotton rolls.

2.2 Precautions

This product is intended to be used only as specifically outlined in these Instructions for Use.

Any use of this product inconsistent with these Instructions for Use is at the discretion and sole responsibility of the dental practitioner.

- Use protective measures for the dental team and patients such as glasses and rubber dam in accordance with local best practice.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the material. Use adequate isolation such as rubber dam.
- Only use the canulas supplied by DENTSPLY for placement of the material.
- Canulas are used for direct intraoral application of the material. Canulas are for single use only, discard in order to prevent cross-contamination. Canula may clog as material may set inside.
- Tightly close syringes immediately after use. Do not allow canula to remain as a syringe cap. Replace original cap tightly after each use.
- The material should extrude easily. **DO NOT USE EXCESSIVE FORCE.** Excessive pressure may cause material and not fully seated canulas to eject from the syringe.
- Interactions:
 - Do not use eugenol- and hydrogen peroxide-containing materials in conjunction with this product since they may interfere with hardening and cause softening of the product.

2.3 Adverse reactions

- Eye contact: Irritation and possible corneal damage.
- Skin contact: Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- Contact with Mucous Membranes: Inflammation (see Warnings).

2.4 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store in sealed packs at temperatures between 10 °C and 24 °C (50 °F and 75 °F).
- Keep out of direct sunlight and protect from moisture.
- Do not freeze.
- Do not use after expiration date.

Humidity can adversely affect the properties of unsealed syringes. For unsealed syringes the following is recommended:

- Store in a dry environment (< 70% relative humidity).
- Use unsealed syringes within 5 months.

3 Step-by-step instructions

3.1 Cleaning

Clean surfaces to be sealed with oil free paste such as Nupro® Prophy Paste. Rinse well with water.

3.2 Isolation

Isolate the teeth to be sealed with rubber dam or cotton rolls and dry each tooth with air, free of oil or water contamination.

3.3 Application of conventional acid conditioner

Refer to adhesive manufacturer's Instructions for Use for tooth conditioning and/or adhesive application.

Dyract® Seal may be used in two different application techniques.

3.3.1 Conventional Technique

1. Application of Phosphoric Acid Conditioner

Apply 35-37% phosphoric acid liquid or gel (e.g. DETREY® Conditioner 36) to the fissures and condition enamel for at least 30 seconds.

2. Rinsing and Drying

Remove acid liquid/gel with aspirator tube and/or vigorous water spray and rinse for at least 15 seconds. Dry surface with a blow of air. Properly etched enamel has a dull, frosty white appearance. If this is not the case, re-etch enamel, rinse and dry as described above.

Once the fissure has been properly treated, it must be kept uncontaminated. If saliva contamination occurs, thoroughly clean with vigorous water spray, dry, and repeat conditioning procedure for 5 seconds only. Rinse and dry as described above.

3. For the complete obliteration of deep and narrow fissures (Total Seal™ Technique) optionally Prime&Bond® NT adhesive may be used (see 3.3.2 Step 2).

3.3.2 Application of NRC™ Non-Rinse Conditioner

Reduced retention time.

The Use of NRC™ Non-Rinse Conditioner instead of conventional acid conditioning results in a reduced lifespan of fissure sealing. Use NRC™ Non-Rinse Conditioner only in the case of poor patient compliance.

1. Application of NRC™ Non-Rinse Conditioner

1.1 Dispense NRC™ Non-Rinse Conditioner into a dappen dish.

1.2 Apply 1 drop of NRC™ Non-Rinse Conditioner with an Applicator Tip or disposable brush to the fissure. Leave undisturbed for 20 seconds. **Do not rinse.**

1.3 Remove excess NRC™ Non-Rinse Conditioner by blowing gently with an air syringe. Once the fissure has been properly treated, it must be kept uncontaminated. If salivary contamination occurs, thoroughly clean with forceful water-spray and repeat the application of NRC™ Non-Rinse Conditioner.

2. Application of Prime&Bond® NT

- 2.1 Dispense Prime&Bond® NT directly onto a fresh Applicator Tip or onto a disposable brush. Alternatively, dispense into a fresh, acetone resistant dappen dish.
- 2.2 Immediately apply 1 drop of Prime&Bond® NT to the fissure.
- 2.3 Leave undisturbed for 20 seconds.
- 2.4 Remove solvent by blowing gently with air from a dental syringe for at least 5 seconds.

3.4 Application

Apply material immediately after enamel pre-treatment.

1. Remove cap from the end of the syringe. To assure free flow of material from syringe, express a small amount onto a mixing pad.
2. Attach disposable steel canula¹ to the end of the syringe. Turn canula clockwise until fully seated. Tug on canula to assure locking into syringe collar.
3. Dispensing tip may be bent for easy access. Gentle curving of the tip is preferable to an acute bend. Bending may occlude the canula. Be sure to check that this has not occurred by extruding a small amount through the tip.
4. Material should flow freely with GENTLE PRESSURE. DO NOT USE EXCESSIVE FORCE. If not, remove syringe from patient field and check for obstruction.
5. Apply material through the canula into the fissure.
6. Pull back syringe plunger slightly after use to prevent excessive flow of material.

Injury due to excessive force.

1. Apply slow and steady pressure on the syringe.
2. Do not use excessive force. Excessive pressure may cause material and not fully seated canulas to eject from the syringe.

3.5 Light curing

1. Light cure each area of sealed surface with a suitable curing light designed to cure materials containing camphorquinone (CQ) initiator, i.e. spectral output containing 470 nm. Minimum light output must be at least 800 mW/cm² exposure for at least 20 seconds. Refer to curing light manufacturer's Instructions for Use for compatibility and curing recommendations.
2. Cure each layer according to the table below.

Inadequate polymerization due to insufficient curing.

1. Check compatibility of curing light.
2. Check curing cycle.
3. Check curing output.

¹ Use the disposable steel canulas included in the Dyract® Seal package. DENTSPLY ReOrder No. of the disposable steel canulas 615.04.707.

Curing Light Output ≥ 800 mW/cm²

Shades	Increment/layer	Curing time
clear	2 mm	10 seconds
white-opaque	2 mm	30 seconds

3.6 Retention and Occlusion Check

1. Remove soft (oxygen-inhibited) surface layer after light curing with cotton pellets or cotton rolls.
2. Check with an explorer for complete coverage and retention.
3. Check with articulating paper for premature occlusal contacts and correct, if necessary, with a finishing bur.

4 Hygiene

4.1 Protecting and disinfecting the syringes

1. Tightly close syringes immediately after use. Do not allow canula to remain as a syringe cap. Replace original cap tightly after each use.
2. To prevent syringes from exposure to spatter or spray of bodily fluids, contaminated hands or oral tissues, use of a protective barrier is recommended.
3. Disinfect contaminated syringes with a water-based hospital-level disinfection solution according to national/local regulations.

Damage to the label by repeated disinfection.

Avoid contamination of syringes.

4.2 Steel canulas – single use

Cross-contamination of steel canulas.

Do not reuse steel canulas.

Properly dispose steel canulas in accordance with local regulations.

5 LOT number () and expiration date ()

1. Do not use after expiration date. ISO standard is used: “YYYY-MM”
2. The following numbers should be quoted in all correspondence:
 - Reorder Number
 - Lot number
 - Expiration date

© DENTSPLY DeTREY 2010-03-01

[These Instructions for Use are based on Master Version 01]

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DENTSPLY DeTREY GmbH
DENTSPLY Limited
DENTSPLY Italia S.r.l.