[ENGLISH] INSTRUCTIONS FOR USE

I. INTRODUCTION

"CLEARFIL Universal Bond Quick" consists of BOND and K-ETCHANT Syringe. The BOND is a light curing bonding agent that allows the treatment of dentin, enamel, and prosthetic materials. Depending on the indication, the BOND is used as self-etching or with K-ETCHANT Syringe for selective enamel etching or total-etching procedures. The BOND is intended to be used for both direct and indirect restorations. "CLEARFIL DC Activator" activates the dual-cure mechanism of the BOND; however, the addition of "CLEARFIL DC Activator" to the adhesive is not required when using it with "CLEARFIL DC CORE PLUS" or "PANAVIA SA Cement Plus". The BOND is available in both Bottle and Unit Dose delivery systems. K-ETCHANT Syringe is an etching gel that consists of 35% phosphoric acid aqueous solution and colloidal silica.

II. INDICATIONS

"CLEARFIL Universal Bond Quick" is indicated for the following uses:

- Direct restorations using light-cure composite resin
- Cavity sealing as a pretreatment for direct restorations
- Treatment of exposed root surfaces
- Treatment of hypersensitive teeth
- Internal repairs of restorations
- Post cementation and core build-ups
- Cementation of indirect restorations

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to methacrylate monomers and this product

IV. POSSIBLE SIDE EFFECTS

[1] The oral mucosal membrane may turn whitish when contacted by the BOND due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.

[2] K-ETCHANT Syringe may cause inflammation or erosion due to its chemistry.

V. INCOMPATIBILITIES

[1] Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.

[2] Do not use histomastics containing ferric compounds, since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva, due to remaining ferric ions.

[3] When using histomastics containing aluminum chloride, minimize the quantity; use caution to prevent contact with the adherent surface. Failure to do so might weaken the bond strength to the tooth structure.

[4] Avoid looking directly at the dental curing light when curing the BOND.

[5] Avoid using the same BOND dispensed into a well of the dispensing dish. Unit Dose and applicator brush for different patients to prevent cross contamination. Unit Dose and applicator brushes are single use only. Discard them after use. The needle is to single-use only. Do not reuse it to prevent cross-contamination. Discard after use.

[6] Wear gloves to take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components in the product.

[7] If the instruments associated with this product are damaged, use caution and protect yourself; immediately discontinue use.

VI. PRECAUTIONS

[1] This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components in the product.

[2] If the patient demonstrates a hypersensitivity reaction such as rash, eczema, features of inflammation, ulcer, swelling, itching, or numbness, discontinue use of the product, remove the product and seek medical attention.

[3] Use caution to prevent the product from coming in contact with the skin or getting into patient's eyes. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.

[4] If the product comes in contact with human body tissues, take the following actions:

-If the product gets in the eyes: Immediately wash the eye with copious amounts of water and consult a physician.

-If the product comes in contact with the skin or the oral mucosa: Immediately wipe the area with a cotton pellet or gauze pad moistened with alcohol and rinse with copious amounts of water.

[5] Use caution to prevent the patient from accidentally swallowing the product.

[6] Avoid looking directly at the dental curing light when curing the BOND.

[7] Avoid using the same BOND dispensed into a well of the dispensing dish. Unit Dose and applicator brush for different patients to prevent cross contamination. Unit Dose and applicator brushes are single use only. Discard them after use. The needle is to single-use only. Do not reuse it to prevent cross-contamination. Discard after use.

[8] Wear gloves to take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components in the product.

[9] If the instruments associated with this product are damaged, use caution and protect yourself; immediately discontinue use.

2. Handling and manipulation precautions

[Common precautions]

[1] This product must not be used for purposes other than specified in the instructions.

[2] The use of this product is restricted to licensed dental professionals.

[3] Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.

[BOND]

[1] The BOND contains ethanol, a flammable substance. Do not use it near an open flame.

[2] To prevent poor performance or poor handling characteristics, observe the specified light-curing times and other handling requirements.

[3] Clean the cavity sufficiently to prevent poor bonding. If the adherent surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.

[4] **Bottle:** Use the light-blocking plate to avoid exposing the material to an operating light or ambient light, and use within the working time indicated below after dispensing or mixing.

Material	Working time
BOND	7 minutes
BOND + CLEARFIL DC Activator	90 seconds

As the ethanol in BOND evaporates, the viscosity increases, thereby making it difficult to apply.

Unit Dose: Apply BOND with the applicator brush immediately after snapping off the container cap.

[5] During the application of BOND to the adherent surface with rubbing motion, move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to the operating light. In addition, after application of BOND, dry sufficiently until BOND does not move by blowing air.

BOND contains ethanol and water. "CLEARFIL DC Activator" contains ethanol. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND or the mixture of BOND and "CLEARFIL DC Activator" does not move. Otherwise, the adhesion effect will be impaired. In order to dry sufficiently, adjust the air pressure according to the shape and size of the cavity and the prosthetic appliance. Use a vacuum aspirator to prevent BOND or the mixture from scattering.

[6] If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat again with BOND.

[7] **Bottle:** Do not mix BOND with other bonding agents except "CLEARFIL DC Activator."

[8] **Unit Dose:** Do not mix BOND with other bonding agents. Do not use Unit Dose with CLEARFIL DC Activator due to the design of the Unit Dose container.

[9] **Unit Dose:** When using for core build-up or cementation, use BOND only with "CLEARFIL DC CORE PLUS" or "PANAVIA SA CEMENT Plus."

[10] **Light-cure:** The mixture of the BOND and "CLEARFIL DC Activator," otherwise, the working time will be shortened dramatically.

[11] **Bottle:** The container should be tightly capped immediately after use to reduce the evaporation of the volatile solvent (ethanol contained in the BOND). If the liquid does not flow easily from the nozzle, do not dispense the BOND forcibly from the plugged container.

[12] **Bottle:** If the BOND has been used with alcohol both before and after use, the BOND may not flow easily, shake the container before use.

[K-ETCHANT Syringe]

[1] Be careful not to contaminate it with saliva or blood. If the treated surface is contaminated, re-treat.

[2] Be careful to avoid cross-contamination. Disinfect the syringe by wiping it with an absorbent cotton with alcohol both before and after use. Cover the entire syringe with a disposable plastic barrier to prevent saliva and blood contamination.

3. If the product adheres to clothing, wash it off with water.

[4] After each use, remove the needle tip from the syringe and recap the syringe immediately and tightly.

[5] Etching the vital dentin may cause post-operative sensitivity.

[Dental light-curing unit]

[1] Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.

[2] The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

3. Storage precautions

[1] The product must be used by the expiration date indicated on the package.

[2] Apply BOND must be refrigerated (2-8°C/36-46°F) when not in use, and should be brought to room temperature for more than 15 minutes before using. Especially, Bottle must be left standing until it comes to room temperature after taking out of the refrigerator; otherwise, an excessive amount of liquid might be dispensed or the liquid may ooze after use.

[3] K-ETCHANT Syringe must be stored at 2 ± 25°C/36 ~ 77°F when not in use.

[4] Keep away from extreme heat, direct sunlight or a flame.

[5] The product must be stored in a proper place where only dental practitioners can access.

VI. COMPONENTS

Please see the outside of the package for contents and quantity.

<Principal ingredients>

[1] BOND

- 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)

- Bisphenol A diglycidylmethacrylate (Bis-GMA)

- 2-Hydroxyethyl methacrylate (HEMA)

- Hydrophilic amide monomers

- Colloidal silica

- Silane coupling agent

- Calcium hydroxide material

- Di-Camphorquinone

- Ethanol

- Water

[2] K-ETCHANT Syringe

- Phosphoric acid

- Colloidal Silica

- Pigment

3. Accessories

- Applicator brush (fine «silver»)

- Dispensing dish

- Light-blocking plate

- Needle tip (E)

III. CLINICAL PROCEDURES

A. Standard procedure I

[1] **Direct restorations using light-cure composite resin**

[2] **Cavity sealing as a pretreatment for indirect restorations**

B. Standard procedure II

[1] **Direct restorations using light-cure composite resin**

[2] **Cavity sealing as a pretreatment for indirect restorations**

C. Treatment of exposed root surfaces

[4] **Treatment of hypersensitive teeth**

-If the product gets in the eyes: Immediately wash the eye with copious amounts of water and consult a physician.

-If the product comes in contact with the skin or the oral mucosa: Immediately wipe the area with a cotton pellet or gauze pad moistened with alcohol and rinse with copious amounts of water.

[5] Use caution to prevent the patient from accidentally swallowing the product.

[6] Avoid looking directly at the dental curing light when curing the BOND.

[7] Avoid using the same BOND dispensed into a well of the dispensing dish. Unit Dose and applicator brush for different patients to prevent cross contamination. Unit Dose and applicator brushes are single use only. Discard them after use. The needle is to single-use only. Do not reuse it to prevent cross-contamination. Discard after use.

[8] Wear gloves to take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components in the product.

[9] If the instruments associated with this product are damaged, use caution and protect yourself; immediately discontinue use.

A-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

A-2. Cavity or root surface preparations

Remove any infected dentin and prepare the cavity in the usual manner. When treating hypersensitive teeth, clean the root surface as usual. Clean thoroughly by spraying with water, and dry with air with cotton pellets.

A-3. Pulp protection

Any actual or near pulp exposure can be covered with a hard setting zinc phosphate cement material. There is no need for cement lining or base. Do not use eugenol materials for pulp protection.

A-4. Pretreatment of tooth

Choose one of the three etching procedures below after applying BOND. [NOTE]

When treating hypersensitive teeth, select section A-4a before applying BOND.

A-4a. Self-etching procedure

Move to section A-5 without etching with K-ETCHANT Syringe.

A-4b. Selective enamel etching procedure

Apply K-ETCHANT Syringe to the uncured and/or cut enamel. Leave it in place for 10 seconds, then rinse and dry.

A-4c. Total-etching procedure

Apply K-ETCHANT Syringe to the entire cavity (enamel and dentin), leave it in place for 10 seconds, then rinse and dry.

A-5. Application of BOND

[1] **Bottle:** Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application. [CAUTION]

Use the light-blocking plate to avoid exposing the material to an operating light or ambient light; use within 7 minutes after dispensing. [NOTE]

Unit Dose: Snap off the container cap.

Use caution not to allow saliva or exudate to contact the treated surface.

[2] Apply BOND with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required. [NOTE]

Use caution not to allow saliva or exudate to contact the treated surface.

[3] Dry the entire cavity wall sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[4] Light-cure BOND with a dental curing unit. (See table "Dental curing unit and curing time".)

Working time will be shortened dramatically when not light-curing the mixture on the adherent surface. [NOTE]

Use caution not to allow saliva or exudate to contact the treated surface.

A-6. Placement of composite resin restorative or treatment of hypersensitive teeth

A-6a. Direct restorations using light-cure composite resin

Place composite resin (e.g. CLEARFIL MAJESTY ES-2, CLEARFIL MAJESTY ES Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

A-6b. Cavity sealing and treatment of exposed root surfaces

Place a thin coat of composite resin (e.g. CLEARFIL MAJESTY ES Plus) flow into the tooth, light-cure according to the manufacturer's instructions. Remove un-polymerized layer with a cotton pellet or a gauze moistened with alcohol.

A-6c. Treatment of hypersensitive teeth

Remove the un-polymerized layer of BOND with a cotton pellet or a gauze moistened with alcohol.

B. Standard procedure II

[3] Intraoral repair of fractured restorations

A. Standard procedure I

Roughen the adherent surfaces using a diamond point or blasting with 30 to 50 µm alumina powder at air pressure of 0.1-0.4 MPa (14-58 PSI / 1-4 kg/cm²). The air pressure should be properly adjusted to suit the material and/or shape of the prosthetic restoration; using caution to prevent it from chipping. Place a bevel at the marginal area.

B-2. K-ETCHANT Syringe treating of adherent surfaces

Apply K-ETCHANT Syringe to the adherent surface (including tooth stripe). Leave in place for 5 seconds, then rinse and dry.

B-3. Application of BOND

Apply BOND with a rubbing motion to the entire adherent surface. Refer to section A-5. [NOTE]

For optimal performance, apply a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER PLUS) to the surface of the silica-based glass ceramic (e.g. conventional porcelain, Lithium disilicate) and apply the silane coupling agent (e.g. ALLOY PRIMER) to the surface including precure metal according to the manufacturer's instructions BEFORE applying BOND.

B-4. Place composite resin restorative

Place composite resin (e.g. CLEARFIL MAJESTY ES-2) into the cavity, light-cure, finish and polish according to the manufacturer's instructions. [NOTE]

Use an opaque resin (e.g. CLEARFIL ST OPAQUE) before composite resin placement to mask metal color.

C. Standard procedure III

[6] Post cementation and core build-up

When using with "CLEARFIL DC CORE PLUS", the use of "CLEARFIL DC Activator" is not necessary. [NOTE]

C-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

C-2. Preparing the root canal

Prepare and clean the root canal in the usual manner.

C-3. Post preparation

Choose either C-3a or C-3b based on the post used. Please follow the instructions for Use of the restorative material. In the absence of specific instructions, we recommend the following procedure:

C-3a. For Glass Fiber Posts

Apply K-ETCHANT Syringe to the entire surface of the post. Leave it in place for 5 seconds, then rinse and dry. [CAUTION]

Do not blast glass fiber posts with alumina powder, or the posts could be damaged.

Avoid all contamination of the surfaces to be treated during pre-treatment and until the final core build-up.

C-3b. For Metal Posts

Roughen the adherent surface by blasting with 30 to 50 µm alumina powder at air pressure of 0.2-0.4 MPa (29-58 PSI / 2-4 kg/cm²). The air pressure should be properly adjusted to suit the material. After blasting, clean the prosthetic restoration by using ultrasonic cleaner for 2 minutes followed by drying it with an air stream.

C-4. Post surface treatment

Choose either procedure based on the material used.

C-4a. When using CLEARFIL DC CORE PLUS

Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.

Apply BOND to the entire post surface with the applicator brush. No waiting time is required.

Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE] To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

C-4b. When using with other dual-/self-cured resin core materials or dual-/self-cured resin cements

[1] **Bottle:** Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them with the applicator brush. [CAUTION]

Use the light-blocking plate to avoid exposing the material to an operating light or ambient light; use within 30 seconds after mixing. **Unit Dose:** Not applicable for this use.

[2] Apply the mixture to the adherent surface.

[3] Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the mixture from scattering.

[NOTE] To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

C-4c. When using with CLEARFIL DC CORE PLUS

Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.

Apply BOND to the entire post surface with the applicator brush. No waiting time is required.

Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent the mixture from scattering.

[NOTE] To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

D-4. Pretreatment of tooth

If necessary, apply K-ETCHANT Syringe. Refer to section A-4.

D-5. Bonding

Choose either procedure based on the material used.

D-5a. When using with PANAVIA SA Cement Plus

Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.

Apply BOND with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required.

Use caution not to allow saliva or exudate to contact the treated surface.

[2] Apply the mixture to the entire post surface with the applicator brush. No waiting time is required.

Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE] To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

D-5b. When using with other dual-/self-cured resin cement

[1] Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them with the applicator brush. [CAUTION]

Use the light-blocking plate to avoid exposing the material to an operating light or ambient light; use within 90 seconds after mixing. **Unit Dose:** Not applicable for this use.

[2] Apply the mixture to the entire cavity wall with the applicator brush. No waiting time is required.

Use caution not to allow saliva or exudate to contact the treated surface.

[3] Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE] To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

D-5c. Pretreatment of tooth

Apply K-ETCHANT Syringe, as needed. Refer to section A-4.

D-6. Bonding

Choose either procedure based on the material used.

D-6a. When using with "CLEARFIL DC CORE PLUS"

[1] Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.

Apply BOND with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required.

Use caution not to allow saliva or exudate to contact the treated surface.

[2] Apply the mixture with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required.

Use caution not to allow saliva or exudate to contact the treated surface.

[3] Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE] To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

D-6b. When using with other dual-/self-cured resin core materials or dual-/self-cured resin cements

[1] Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them with the applicator brush. [CAUTION]

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