



NEW

radii plus + attachments

your choice of
accessories

SDI

full arch
bleach LED
attachment

for in-office
bleaching of full arch

single tooth
bleach LED
attachment

for in-office
bleaching of
single tooth

diagnostic
LED
attachment

see your patients
in a different light!

orthodontic
curing LED
attachment

fast and efficient
bonding of brackets

contra angle
curing LED
attachment

a decisive step
forward to your
curing needs

right angle
curing LED
attachment

great for
occlusal curing

Please note: Your new Radii Plus "+" LED curing light will work only with the new Radii Plus "+" attachments.

Look for the small "+" sign on the Control section which will indicate the improved Radii Plus.

full arch bleach LED attachment

for in-office bleaching

- Beeps every minute
- Heats to 37°C / 98.76°F



order: 5600092

bleach arch stand

lightweight & compact

- fits easily into a cupboard



order: 5600095

single tooth bleach LED attachment

for in-office bleaching

- Beeps every 30 seconds
- Heats to 45°C / 113°F



NEW

order: 5600258

contra angle LED attachment

a decisive step forward
to your curing needs

- Just like your handpiece
- Permits access to areas difficult to cure



NEW

order: 5600276

right angle LED attachment

curing tip

- Ease of access
- Great for occlusal curing



NEW

order: 5600260

orthodontic LED attachment

fast and efficient
bonding of brackets

- Beeps every 10 seconds
- Runs continuously for 5 minutes



NEW

order: 5600259

diagnostic LED attachment

examination using
transillumination

- Lights up those hard to see areas to assist in identification of cracks & decay



NEW

order: 5600268

LED radiometer

instant digital intensity
reading

- specially calibrated for any LED curing light



order: 5600028

curing tacking tip

great for seating of
veneers

- place over the lens cap when curing
- 1mm and 2mm openings
- Autoclavable



order: 5600079

diagnostic adapter tip

produces a fine beam of
transilluminating light

- replace lens cap with this adapter tip when using the diagnostic attachment
- Autoclavable



order: 5600078